



Title: UV Light Exposure

Principle:

To prevent or reduce ultraviolet light exposure from instrument sources for all personnel working in Anatomic Pathology. UV light may cause corneal or skin burns from direct or deflected light sources

Procedure:

In Anatomic Pathology the only source of UV light is Biological Safety Cabinets. Signs will be posted on Biological Safety Cabinets that this device produces potentially harmful ultraviolet (UV) light. Protect eyes and skin from exposure while light is switched on.

Written by: Gail Zander, CT (ASCP), 8/25/2014

Approval of Procedure:

Medical Director Signature: *m. Matt Zucchi, MD.*

Date: *8/26/14*

